

**Bedsore Free Device
Order Form**
<http://www.bedsorefree.net/>

| | | |
|-------|----------|--|
| _____ | \$496.50 | BFD-A1 - 35" Bedsore Free Device comes with 30 ft. of 3/8" I.D. vinyl tubing and hose plugs. (Needs a low pressure air/vacuum pump to work- See BFD-J1) (The Device can be secured by standard cord or nylon string or by purchasing one of our support straps.(BDF-W1 or BFD-W2)) |
| _____ | \$546.50 | BFD-A2 - 42" Bedsore Free Device comes with 30 ft. of 3/8" I.D. vinyl tubing and hose plugs. (Needs a low pressure air/vacuum pump to work- See BFD-J1) (The Device can be secured by standard cord or nylon string or by purchasing one of our support straps.(BDF-W1 or BFD-W2)) |
| _____ | \$596.50 | BFD-A3 - 48" Bedsore Free Device comes with 30 ft. of 3/8" I.D. vinyl tubing and hose plugs. (Needs a low pressure air/vacuum pump to work- See (BFD-J1) (The Device can be secured by standard cord or nylon string or by purchasing one of our support straps.(BDF-W1 or BFD-W2)) |
| _____ | \$536.50 | BFD-B1 - 35" Bedsore Free Device comes with 30 ft. of 3/8" I.D. vinyl tubing, hose plugs and Electric Air/Vacuum Pump.(BFD-J1) (The Device can be secured by standard cord or nylon string or by purchasing one of our support straps.(BDF-W1 or BFD-W2)) |
| _____ | \$586.50 | BFD-B2 - 42" Bedsore Free Device comes with 30 ft. of 3/8" I.D. vinyl tubing, hose plugs and Electric Air/Vacuum Pump.(BFD-J1) (The Device can be secured by standard cord or nylon string or by purchasing one of our support straps.(BDF-W1 or BFD-W2)) |
| _____ | \$636.50 | BFD-B3 - 48" Bedsore Free Device comes with 30 ft. of 3/8" I.D. vinyl tubing, hose plugs and Electric Air/Vacuum Pump.(BFD-J1) (The Device can be secured by standard cord or nylon string or by purchasing one of our support straps.(BDF-W1 or BFD-W2)) |

| | | |
|-------|----------|---|
| _____ | \$645.50 | BFD-C1 - 35" Bedsore Free Device comes with 30 ft. of 3/8" I.D. vinyl tubing and with the Manual 10-Valve Control Unit (BFD-L1) (Needs non-hand-held low pressure air/vacuum pump to work - See BFD-K1) (The Device can be secured by standard cord or nylon string or by purchasing one of our support straps.(BDF-W1 or BFD-W2)) |
| _____ | \$695.50 | BFD-C2 - 42" Bedsore Free Device comes with 30 ft. of 3/8" I.D. vinyl tubing and with the Manual 10-Valve Control Unit (BFD-L1) (Needs non-hand-held low pressure air/vacuum pump to work - See BFD-K1) (The Device can be secured by standard cord or nylon string or by purchasing one of our support straps.(BDF-W1 or BFD-W2)) |
| _____ | \$745.50 | BFD-C3 - 48" Bedsore Free Device comes with 30 ft. of 3/8" I.D. vinyl tubing and with the Manual 10-Valve Control Unit (BFD-L1) (Needs non-hand-held low pressure air/vacuum pump to work - See BFD-K1) (The Device can be secured by standard cord or nylon string or by purchasing one of our support straps.(BDF-W1 or BFD-W2)) |
| _____ | \$846.50 | BFD-D1 - 35" Bedsore Free Device comes with 30 ft. of 3/8" I.D. vinyl tubing and with the Manual 10-Valve Control Unit (BFD-L1) and with Air/Vacuum Pump (BFD-K1). (The Device can be secured by standard cord or nylon string or by purchasing one of our support straps.(BDF-W1 or BFD-W2)) |
| _____ | \$896.50 | BFD-D2 - 42" Bedsore Free Device comes with 30 ft. of 3/8" I.D. vinyl tubing and with the Manual 10-Valve Control Unit (BFD-L1) and with Air/Vacuum Pump (BFD-K1). (The Device can be secured by standard cord or nylon string or by purchasing one of our support straps.(BDF-W1 or BFD-W2)) |
| _____ | \$946.50 | BFD-D3 - 48" Bedsore Free Device comes with 30 ft. of 3/8" I.D. vinyl tubing and with the Manual 10-Valve Control Unit (BFD-L1) and with Air/Vacuum Pump (BFD-K1). (The Device can be secured by standard cord or nylon string or by purchasing one of our support straps.(BDF-W1 or BFD-W2)) |

| | | |
|-------|-----------|---|
| _____ | \$1175.50 | BFD-E1 - 35" Bedsore Free Device comes with 30 ft. of 3/8" I.D. vinyl tubing, the Manual Switch 10-Solenoid Valve Control Unit (BFD-N1) and the Air/Vacuum Pump (BFD-K1) (The Device can be secured by standard cord or nylon string or by purchasing one of our support straps.(BDF-W1 or BFD-W2)) |
| _____ | \$1215.50 | BFD-E2 - 42" Bedsore Free Device comes with 30 ft. of 3/8" I.D. vinyl tubing, the Manual Switch 10-Solenoid Valve Control Unit (BFD-N1) and the Air/Vacuum Pump (BFD-K1) (The Device can be secured by standard cord or nylon string or by purchasing one of our support straps.(BDF-W1 or BFD-W2)) |
| _____ | \$1255.50 | BFD-E3 - 48" Bedsore Free Device comes with 30 ft. of 3/8" I.D. vinyl tubing, the Manual Switch 10-Solenoid Valve Control Unit (BFD-N1) and the Air/Vacuum Pump (BFD-K1) (The Device can be secured by standard cord or nylon string or by purchasing one of our support straps.(BDF-W1 or BFD-W2)) |
| _____ | \$1725.50 | BFD-F1 - 35" Bedsore Free Device comes with 30 ft. of 3/8" I.D. vinyl tubing, with the Manual Switch 10-Solenoid Valve Control Unit (BF-N1), with the Automatic Add On Unit (BFD-P1) and with the Air/Vacuum Pump (BFD-K1), includes Software to control device automatically. (Requires a computer with 2 available comports) (The Device can be secured by standard cord or nylon string or by purchasing one of our support straps. (BDF-W1 or BFD-W2)) |
| _____ | \$1765.50 | BFD-F2 - 42" Bedsore Free Device comes with 30 ft. of 3/8" I.D. vinyl tubing, with the Manual Switch 10-Solenoid Valve Control Unit (BF-N1), with the Automatic Add On Unit (BFD-P1) and with the Air/Vacuum Pump (BFD-K1), includes Software to control device automatically. (Requires a computer with 2 available comports) (The Device can be secured by standard cord or nylon string or by purchasing one of our support straps. (BDF-W1 or BFD-W2)) |
| _____ | \$1799.50 | BFD-F3 - 48" Bedsore Free Device comes with 30 ft. of 3/8" I.D. vinyl tubing, with the Manual Switch 10-Solenoid Valve Control Unit (BF-N1), with the Automatic Add On Unit (BFD-P1) and with the Air/Vacuum Pump (BFD-K1), includes Software to control device automatically. (Requires a computer with 2 available comports) (The Device can be secured by standard cord or nylon string or by purchasing one of our support straps. (BDF-W1 or BFD-W2)) |

| | | |
|-------|-----------|--|
| _____ | \$45.00 | BFD-J1 - 120VAC Electric Air/Vacuum Pump - made to hold in your hand while you air up or deflate inflatable devices. Noisy while running. |
| _____ | \$245.00 | BFD-K1 - Electric Air/Vacuum Pump-Stand Alone (low decibels) |
| _____ | \$160.00 | BFD-L1 - Manual 10-Valve Control Unit without Air/Vacuum Pump. This Unit needs to sit on a support unit in front of the foot of the bed at a height that would be comfortable to the care provider. It can be set on the floor and slid under the bed when not needed. |
| _____ | \$389.00 | BFD-M1 - Manual 10-Valve Control Unit with the Air/Vacuum Pump (BFD-K1). This Unit needs to sit on a support unit in front of the foot of the bed at a height that would be comfortable to the care provider. It can be set on the floor and slid under the bed when not needed. |
| _____ | \$485.00 | BFD-N1 – Manual Switch 10-Solenoid Valve Control Unit without Air/Vacuum Pump. Designed to slip under most beds out of the way. The switch box used to control the Bedsore Free Device can be used by the patient or care provider. (Needs a external air and vacuum supply to operate unit.) |
| _____ | \$699.00 | BFD-O1 – Manual Switch 10-Solenoid Valve Control Unit with the Air/Vacuum Pump (BFD-K1). The switch box used to control the Bedsore Free Device can be used by the patient or care provider. |
| _____ | \$588.50 | BFD-P1 – Automatic Add On Unit for BFD-N1 & BFD-O1, includes software to control the device automatically. (Requires computer with 2 available comports) This unit contains relay boards used to energize and de-energize the electric solenoids found in the control unit. |
| _____ | \$1250.00 | BFD-S1 - Device comes with the Manual Switch 10-Solenoid Valve Control Unit (BFD-N1), with the Automatic Relay Unit (BFD-P1) and with the Air/Vacuum Pump (BDF-K1), includes software to control the device automatically. (Requires computer with 2 available comports) |

